



ZIRCONIA PER STRUTTURE PROTESICHE DENTALI

Composizione chimica		DB+™	DB+LUX™
• Y ₂ O ₂	wt%	5,2	5,2
• Al ₂ O ₃	wt%	0,25	0,05
• SiO ₂	wt%	<0,02	<0,02
• Fe ₂ O ₃	wt%	<0,01	<0,01
• Particle Size	nm	90	40
• Crystallite Size	nm	36	36
• LOI	wt%	3,4	4,0
• BET	m ² /g	7	13
• Binder	wt%	3	3
• Bulk Density	g/cm ³	1,2	1,2

BETTINI S.p.A.

Via Industriale, 11 – 23804 MONTE MARENZO (LC) * TEL : 0341 – 631.588 FAX : 0341 – 631.589

* P. Iva e C.F. 08212540150

Società soggetta a direzione e coordinamento da parte di Fintex Srl – Società Unipersonale